

BarbLock Retainers
Fitting Matrix
Reference

Decimal	Size Fitting	A Ø	B Ø	C Ø	D Ø	E Ø	F Ø	G Reference	H	I	J
0.125	1/8"	0.19	0.125	0.09	0.11	0.25	0.15	0.485	0.10	0.175	0.175
0.187	3/16"	0.240	0.187	0.140	0.170	0.350	0.193	0.58	0.10	0.29	0.29
0.250	1/4"	0.310	0.250	0.183	0.213	0.450	0.256	0.70	0.10	0.35	0.35
0.375	3/8"	0.450	0.375	0.285	0.315	0.625	0.382	0.86	0.10	0.43	0.43
0.500	1/2"	0.600	0.500	0.410	0.440	0.750	0.510	1.06	0.10	0.52	0.54
0.750	3/4"	0.928	0.750	0.650	0.680	1.000	0.765	1.48	0.10	0.86	0.62

Equivalent Tubing ID Barb OD Fitting Top Fitting ID Base ID Tube Stop Neck OD Retain OD Overall Neck Hgt. Barb Hgt.

Note	DESCRIPTION	Dimension	
Inside Corners	Minimum .010 Rad.	All	

REV.	DATE	BY	DESCRIPTION
0	12-03-07	AAW	Release for Production

Uncontrolled Document

TOLERANCES UNLESS OTHERWISE NOTED:
TOLERANCES ARE NOT ADDITIVE.

See Matrix For Dimensions

CONFIDENTIAL: NOT TO BE REPRODUCED OR USED IN ANY WAY WITHOUT WRITTEN CONSENT FROM Twin Bay Medical, Inc.
11550 US-31 South
Williamsburg, MICHIGAN

Twin Bay Medical, Inc.
11550 US-31 South
Williamsburg, Michigan 49690
Telephone (231) 254-0101
Fax (231) 254-0100

DESIGNED BY: AAW	DATE: 12-29-07	CHECKED BY: JAW	DATE: 11-29-07	CUSTOMER: BarbLock	PROJECT:
APPROVED BY: AAW	DATE: 12-29-07	NEXT ASSEMBLY:		SCALE: 100% (SEE FITTING)	JOB NO:
NAME: Master Fitting Matrix			DWG NO:	REV: 0	