

**BarbLock Retainers**  
Fitting Matrix  
Reference

Decimal	Size Fitting	A ∅	B ∅	C ∅	D ∅	E ∅	F ∅	G Reference	H	I	J
0.125	1/8"	0.19	0.125	0.09	0.11	0.25	0.15	0.35	0.10	0.175	0.175
0.187	3/16"	0.240	0.187	0.140	0.170	0.350	0.193	0.58	0.10	0.29	0.29
0.250	1/4"	0.310	0.250	0.183	0.213	0.450	0.256	0.70	0.10	0.35	0.35
0.375	3/8"	0.450	0.375	0.285	0.315	0.625	0.382	0.86	0.10	0.43	0.43
0.500	1/2"	0.600	0.500	0.410	0.440	0.750	0.510	1.06	0.10	0.52	0.54
0.750	3/4"	0.928	0.750	0.650	0.680	1.000	0.765	1.48	0.10	0.86	0.62
1.000	1"	1.210	1.020	0.900	0.930	1.460	1.090	1.49	0.10	0.86	0.63
Equivalent	Tubing ID	Barb OD	Fitting Top	Fitting ID	Base ID	Tube Stop	Neck OD	Fitting Hgt.	Foot Thickness	Neck Hgt.	Barb Hgt.

<b>Note</b>	<b>DESCRIPTION</b>	<b>Dimension</b>	
Inside Corners	Minimum .010 Rad.	All	

  

REV.	DATE	BY	DESCRIPTION
0	12-03-07	AAW	Release for Production

Uncontrolled Document

CONFIDENTIAL: NOT TO BE REPRODUCED OR USED IN ANY WAY WITHOUT WRITTEN CONSENT FROM Twin Bay Medical, Inc. 11590 US-31 South Williamsburg, MICHIGAN

TOLERANCES UNLESS OTHERWISE NOTED: TOLERANCES ARE NOT ADDITIVE.

See Matrix For Dimensions

DRAWN BY: AAW DATE: 12-28-07 CHD BY: TAY DATE: 11-28-07  
 APP'D BY: DMW DATE: 12-28-07 MEX AS: 11-28-10

NAME: Master Fitting Matrix

Twin Bay Medical, Inc.

11590 US-31 South  
Williamsburg, Michigan 49690  
Telephone (231) 264-0101  
Fax (231) 264-0100